

## Bakterie-Virus/HME-filter



Bakterie/Virus/HME-filter



Bakterie/virusfilter

### Produktinformation

**ALTERA**<sup>®</sup>

- Komplet sortiment av bakterie-virus/HME-filter.
- Bakterieeffektivitet: 99.9999%, Viruseffektivitet 99.999%.

Material: SBC och PP



### Orderinformation

Artnr	Beskrivning	Information	Avd förp
AL-08016.V001	Bakterie-virus/HME filter vx.	Tidal volym: 150-1500ml Dead space: 53ml Bakterie: 99.9999% Virus: 99.999% Befuktning: 36,8 mgH <sub>2</sub> O/L Vikt: 30gram	50st
AL-08018.V001	Bakterie-virus filter	Tidal volym: 150-1500ml Dead space: 33ml Bakterie: 99.9999% Virus: 99.999% Vikt: 19gram	50st